
Manual of HDAC6 DETECTOR

open-source software


Step 1. Select QSAR models for prediction.

You can choose a model developed using the Morgan fingerprints and gradient boosting method (GBM) or a model developed using topological fingerprints, and multilayer perceptron (MLP) classifier.



Step 2. Select input molecular files.

If you choose smiles, please, directly paste the SMILES representation of the desired chemical structure and press Ctrl+Enter. If the entered chemical structure is correct, the application will generate a 2D image of the studied compound. Click on the "PREDICT COMPOUND FROM SMILES" button.



HDAC6 DETECTOR


The application provides an alternative method for




If the entered structure is incorrect, the application reports an error that has occurred, the details of which can be found through the application manager.

```
rdkit.Chem.rdchem.AtomValenceException: This app has encountered an error. The original error
message is redacted to prevent data leaks. Full error details have been recorded in the logs (if you're
on Streamlit Cloud, click on 'Manage app' in the lower right of your app).
  File "/home/appuser/venv/lib/python3.9/site-packages/streamlit/runtime/scriptr
      exec(code, module.__dict__)
  File "/app/hdac6_detector/HDAC6_detector_app.py", line 182, in <module>
      smiles=standardize_smiles(compound_smiles)
  File "/home/appuser/venv/lib/python3.9/site-packages/molvs/standardize.py", li
      mol = Standardizer().standardize(mol)
  File "/home/appuser/venv/lib/python3.9/site-packages/molvs/standardize.py", li
      Chem.SanitizeMol(mol)
```

To clarify the details of the error, please click here_




If you choose a file *sdf, that may contain a different number of chemical structures, please specify the path to this file on your computer's hard drive. In this case, you need to click the "Browse files" button.



HDAC6 DETECTOR

The application provides an alternative method for assessing the potential of chemicals to be Histone




Step 3. *Prediction results*. The form of presentation of the results depends on the type of descriptors selected, as well as the format of the input chemical data. For example, when selecting SMILES, the results will be displayed for a single molecule. The chemical fragments are colored in green (predicted to reduce inhibitory activity) or magenta (predicted to increase activity HDAC6 inhibitors). The gray isolines separate positive and negative contributions.

Prediction results:


HDAC1: Active

Applicability domain (AD): Inside AD

Predicted fragments contribution:



The Bioavailability Radar: compliance the Lipinski's rule of five




3D structure of the studied compound:


4.0000

10.0000

Number of hydrogen bond acceptors(HBAs)



If you select a file *sdf, the results will be presented in tabular form. If incorrect structures are detected in the file *sdf, the corresponding information will appear in the section "CHEMICAL STRUCTURE VALIDATION AND STANDARDIZATION"



2. RESULTS OF PREDICTION:



The prediction results for correct chemical structures are displayed in a table that can be downloaded, or separately for each molecule if you click the "Show results and map of fragments contribution for each molecule separately" button.

